








Diabact® UBT 50 mg Helicobacter-Test

Ticks all the Boxes

- 
Faster Test
 Samples are ready for analysis in only 10 minutes.
- 
Safe
 Uses the naturally occurring, stable non-radioactive isotope ¹³C.
 Can be used during pregnancy and lactation.
 Approved in Australia by the TGA.
- 
Reliable
 98 - 100% sensitivity and specificity.
- 
More Convenient
 Easy to handle, easy to swallow. Patient takes just a single tablet with water. No pre-test orange juice or citric acid required.
- 
Cost Effective
 Diabact® UBT is significantly less expensive than older tests. Discounts are available for larger orders.



MDA

Diabact® UBT 50mg

The **Fastest**
 Non-Radioactive
UREA BREATH TEST
 to Detect
HELICOBACTER PYLORI

	ENDOSCOPIC	SEROLOGY	LIQUID UBTS	DIABACT® UBT
PATIENT COMFORT AND CONVENIENCE	•	••	••	•••
TIME	•	•••	•	•••
PREPARATION	•	••	••	•••
COST	•	•••	••	•••
SENSITIVITY AND SPECIFICITY	•••	•	••	•••



Medical Diagnostics Australasia

MDA

Medical Diagnostics Australasia
 PO Box 5320
 Clayton Vic 3168

Tel: 03 9551 5105
 e-mail: mdaus@bigpond.com.au
www.medicaldiagnostics.net.au

TGA approved: AUST R 146728 AUST R 135657

Diabact® UBT 50mg

The **Fastest**
 Non-Radioactive
UREA BREATH TEST
 to Detect
HELICOBACTER PYLORI



10 minute test for *Helicobacter*



MDA

Patient Samples in 3 Easy Steps



1 Step 1.

Take the baseline breath samples

2 Step 2.

Take the Diabact® UBT tablet with a glass of water

3 Step 3.

Take the 10 minute breath samples

Clinical Trial Proven Reliability

Clinical trials have demonstrated excellent performance by Diabact® UBT in European and Chinese populations respectively.

Gatta et al 2003.

A rapid, low-dose, ¹³C-urea tablet for the detection of *Helicobacter pylori* infection before and after treatment. (Aliment. Pharmacol. Ther. 17: 793-798)

Wong et al 2003.

A rapid-release 50mg tablet-based ¹³C-urea breath test for the diagnosis of *Helicobacter pylori* infection. (Aliment. Pharmacol. Ther. 17: 253-257)

They conclude: 'New 10 minute ¹³C-UBTs using tablet formulations of ¹³C-urea with citric acid are reliable for the assessment of *H. pylori* status pre- and post-treatment'.

Unique Tablet Formulation

Diabact® UBT 50mg is a TGA approved product manufactured by Kibion AB, a Swedish company specializing in targeted drug delivery systems.

The rapid release tablet is uniquely formulated to prevent oropharyngeal urease contamination of the test and to provide a high concentration of the substrate at the site of the infection.

Diabact® UBT's solid dosage system minimizes false positive results in early breath samples by eliminating urea hydrolysis in the mouth and on the way to the stomach. In addition, the acid supplementation in the tablet reduces the risk of false negative results.

